






### Related Videos

-  Sublimation lyophilic dryer for...  
28.700K views
-  Automatic machine for...  
1.555K views
-  Filling and sealing equipment...  
4.520K views
-  Automatic machine for fillin...  
2.020K views
-  Powder mixing equipment in dru...  
2.171K views
-  Automatic equipment for...  
100.545K views
-  Automatic labeling...  
74.710K views

## High speed automatic powder filling for hard gelatin capsules

1.754K views · 09.09.2020

1.399K likes 50 comments

Screenshots



## Pharmaceutical Tablet Press

597 Videos ·



High speed automatic powder filling for hard gelatin capsules Pharma-Manager.com Online video of production and laboratory equipment from leading manufacturers and suppliers of pharmaceutical industry.

OUR KEYWORDS: Equipment, catalog, pharmaceutical, technology, manufacturing, laboratory, GMP, news, pharmaceuticals, products, industry, plant, design, clean rooms, gelatine capsule, vial, ampoule, tube, cap, capping, labeler, labeller, machine, tablet, press, tabletpress, pharma, filling, dispenser, pharma, packing, blister, marker, labeling, printer, application, manual, capsule filling, capsule filler, sale, buy, photo, characteristics, sell, soldering, sealer, powder, granules, liquid, solid, soft, dry.

Our visitors will be able to see all processes how it is done. We also help you choose different kinds of pharmaceutical equipment – laboratory tablet press, rotary tablet press, hydraulic press, granulators, mixers, drying, vacuum powder conveyor, blistering machines, packaging machines, capsule filling machines and many other kinds of equipment. Our consultants will help with the discussion, selection, payment, and delivery of pharmaceutical equipment to your city. For over 17 years we have been working in this industry for our customers.

Company CEO Mr. Roman Tsibulsky

Write: [info@pharma-manager.com](mailto:info@pharma-manager.com)

Mobile phone WhatsApp +375-29-393-23-82

Skype: RomanTsibulsky

Capsulator type of pharmaceutical equipment. Used to automatically fill gelatin capsules with powders, or to form capsules of gelatinous materials, followed by filling with oil or fat. Capsules can be referred to as manual semi-automatic machines, as well as automatic equipment of high productivity. Capsule production the technological process of obtaining pharmaceutical products in the production of gelatin capsules. The term can be attributed to the process of obtaining hard gelatin capsules, as well as soft gelatin capsules (suture and seamless). Capsules type of pharmaceutical products. Capsules are formed using gelatin, which forms the walls of the shell. They are used for dosage of medicinal preparation in solid and liquid forms. Capsules are water-soluble. Hard gelatin capsule type of pharmaceutical products. Stable thin-walled gelatin bowls made of two parts, bases and lids. They are used for packaging grains and powdered medicines. Solid capsules have standard sizes and are subdivided into numbers 000,00,0,1,2,3,4.

OUR CATALOG OF PHARMACEUTICAL EQUIPMENT: Equipment for preparation and transportation of tablet mass; Early boring and product feeding equipment; Centrifuge Sieve; Vibrating sieve for sieving powders and granules; Cone Pharmaceutical Mills; Pharmaceutical hammer mills; Pellet mixers for industrial production; Fluid Bed Pellet Dryers; Bins for storing tablets, powders and pellets; Bin mixers for powders and granules; Laboratory pharmaceutical equipment; Industrial equipment for medicine production; Lifting columns for containers; Bottle and vial washing; Semi-automatic bottle and vial washer; Automatic Bottle Washing; Tablet presses for tablet pressing; Laboratory rotary tablet press; Laboratory rotary tablet press; Compact rotary tablet press; Washable rotary tablet press; Washable rotary tablet press; Rotary tablet press for powder pressing; Press tool for tablet equipment; Coater for the application of the shell on the tablets; Container washing and sterilization equipment; Semi-automatic bottle washing machine for syringe washing; Automatic washing and drying machine for ampoules, vials and cartridges; A pyrogenic tunnel for sterilization of ampoules, vials, cartridges; Powder filling and dispensing lines; Liquid solution bottling and powder dispensing machines; Bottling machine for products with low viscosity; Bottling machine and double capping bottles; Sunset Machine for Vials of Different Forms; Bottling and capping machine for ophthalmic products; Glass bottle filling and capping machine; Plastic Bottle Closing Machine; Bottling machine in glass ampoules and penicillin bottles; Filling machines in plastic syringes and cartridges; Bottling machine in vials, syringes and cartridges; Syringe and cartridge filling machine; Filling in plastic bags and bottles; Infusion Packaging Machine; Bottling machine for ready-made infusion bags; Line for forming containers from granules, with subsequent filling and capping; Filling and capping machine for PE and PP bottles; Dispersion into plastic bottles and vials; Inspection equipment for drug testing; Inspection machine for glass ampoules and penicillin bottles; Inspection machine for syringes and cartridges; Inspection machine for



infusion packs; Inspection machine for large glass and plastic bottles; Carton machines for blisters and tubes packaging; Automatic drug cartridge machine; Vertical cardboard machine; Horizontal cardboard machine; Folding machine.



Show more

#Adjustment #Budgetary #Cheap #Chile #Comments #Czech Republic #Ecuador #Engineer #Know-how #Mauritania #Norway #Now #On credit #Order #Papua New Guinea #Pay #Pharma #Purchase #Rate #Recommendation #Technology #Trinidad and Tobago #Vietnam



## Products



**Preparation of Gelatin**

36 videos



**Powders**

59 videos



**Powder Mixer**

65 videos



**Powder in bags**

42 videos



**Powder compactor**

39 videos



LAST

High speed automatic rotary tablet press for tablet production

Pharmaceutical Tablet Press • 09.09.2020

Homogenizer to prepare a cream and ointment for cosmetics and pharmaceutical products

NEXT



Pharmaceutical Tablet Press • 09.09.2020

More links



- Vacuum packaging machines and
- Automatic celloperator envelope
- Packaging equipment china
- Dispenser aggressive liquids
- Production lodgements blisters

- Freeze dried drying of vegetables
- Packing list granulator
- Plastic tube price
- Packing medicines in the capsules
- Dietary supplements requirements to the package

- Buy vacuum packaging machine for the house
- Equipment sugar packaging in a stick
- Sizes solid gelatin capsules
- Drawing tablet press
- Molded blister



## You may also like this



Automatic equipment for feeding plastic caps with bottle closure system

Pharmaceutical Packaging • 24.09.2020 • 15.475K views



Automatic equipment for tablet packing in soft aluminum packaging

Pharmaceutical Tablet Press • 23.11.2020 • 45.410K views



Group Cardboard boxes cellophane wrap equipment USA

Pharmaceutical Tablet Press • 04.09.2020 • 1.797K views

## Leave us a comment

Enter Your Name\*

Enter Your Email\*


Enter Your Comment





Post Comment

# PHARMA VIDEO



Demonstration of equipment for the production of tablets, solid gelatin capsules, soft gel capsules, dragee, granules, powder, cream, ointment, suspension, liquid medicines, suppositories, gel.

Los Angeles, CA, U.S.

WhatsApp +375-29-393-23-82

info@pharma-manager.com

## Information

[About Us](#)

[Contact Us](#)

[Catalog](#)

[Report videos](#)

[Popular Videos](#)

[Top 10 photos](#)

Select Country ▾

## Categories

[Blisters](#)

[Tablets](#)

[Tube](#)

[Powder Press](#)

[Powder Mixing](#)

[Vibrating Screen](#)

[Sublimation](#)

## Quick Links

[Soft medical Packing](#)

[Vibrating Screen](#)

[Agar Oil](#)

[Granules](#)

[Machine Sealing](#)

[Equipment for Cleaning](#)

[Vacuum Dispensers](#)

## Latest Video



**Automatic equipment for...**  
10.12.2020 ·  
110.570K views



**Automatic equipment for...**  
10.12.2020 ·  
145.400K views



**Automatic equipment for...**  
10.12.2020 ·  
150.770K views

